de

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BegoPal® 300

Dental alloy

en

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The enclosed BEGO-GOLD card is part of these instructions for use!

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ISO 9693 / ISO 22674

Alloy characteristics

Bio-Certificate	<input checked="" type="checkbox"/>	Melting interval [°C]	1320–1175
Type (ISO 22674)	4	Casting temperature [°C]	1390
BEGO-GOLD colour code	white 8	Preheating temperature [°C]	850
Density [g/cm ³]	11.0	CTE 25–500 °C [10 ⁻⁶ K ⁻¹]	13.8
Average grain size [µm]	35	CTE 20–600 °C [10 ⁻⁶ K ⁻¹]	14.0
Vickers hardness (HV 5)	* 220/240/260	Heat treatment	500 °C 15 min.
Elongation limit (R _{p0,2}) [MPa]	* 450/520/540	Soft annealing	750 °C 10 min.
Ductile yield (A ₅) [%]	* 30/28/25	(then quenching in water at 20 °C)	
Modulus of elasticity [GPa]	approx. 135		

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Richtanalyse in Masse-% (Elemente)

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Safety hint

Metal dust is harmful to your health.
When deflasking and blasting use a suction
extraction system and breathing mask
type FFP3-EN149:2001!

Standard analysis, % by weight (elements)

Pd	75.4
In	6.3
Ag	6.2
Au	6
Ga	6
Ru	

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
**Instructions for use****Modelling:**

- Minimum metal thickness (after grinding): for ceramic veneering 0.4 mm, for acrylic veneering with retention pearls 0.3 mm.
- Connecting parts between the pontics shall be as thick and high as possible (at least 3.5 mm height and 2.5 mm width).
- Provide palatal interdental garlands as reinforcement for wide span bridges.

Investing: Use phosphate bound crown- and bridge investment (e. g. Bellavest® SH, BellaStar XL).

Casting and finishing: General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only.

- In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material.
- Use ceramic crucibles and sprinkle some Auromelt HF melting powder on the casting ingots.
- Keep on heating after complete collapsing of the casting pieces!

Flame melting:	10–15 seconds
HF-induction melting:	approx. 10 seconds (depending on induction capacity)
Resistance heating:	3–4 minutes

Use fine carbide, ceramically bonded stones or BEGO sintered diamond milling tools for finishing.

Ceramic: Use veneering ceramics in accordance with ISO 9693 with firing temperatures of up to approx. 980 °C (e. g. VM 13, Duceram KISS, Omega 900, BeCe Press). Always follow the ceramic manufacturer's instructions!

- Always blast the surface to be veneered (Korox® 110, 2 – max. 3 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.).
- Allow to cool down normally after firing.

Oxide firing:

- Maintain oxide firing at 960 °C without vacuum for 2–3 minutes (Omega 900: 900 °C).
- The oxide should be blasted again prior to application of the ceramics.

Ceramic firing: In the case of washbrand, apply a very thin layer of the paste-like basic substance using a flat brush and allow to predry for 5 minutes. Opaque firing on entire surface.

Soldering:

- Support object in a soldering block of Bellatherm®. Prepare a gap of max. 0.2 mm with parallel walls.
- Soldering before firing with the flame (1125 °C): BegoStar®-Solder (REF 61081) and Minoxyd or Fluxsol.
- Soldering after firing in furnace (810 °C): BEGO-Gold-Solder I (REF 61017) and Minoxyd. Allow to cool normally.
- Acid treatment for residues of Minoxyd in Aurocid (60 °C, approx. 1 minute). Then clean thoroughly (steam clean or boil in aqua dest.).

Laser welding: Filler material: BegoPal® 300 wire diameter 0.4 mm (REF 61165) or accurately fitting, custom-ground casting pieces.

Secondary effects: Like e. g. allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

Reciprocal actions: In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known incompatibilities and allergies to contents of the alloy.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products undergo constant further development and are therefore subject to modification regarding design and composition.

BEGO