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Alloy characteristics

Biocertificate	<input checked="" type="checkbox"/>	Melting interval [°C]	1280–1230
Type	extra-hard (4)	Casting temperature [°C]	1380
BEGO-GOLD colour code	white 8	Preheating temperature [°C]	850
Density [g/cm ³]	13,8	CTE 20–600 °C (25–500 °C)	
Average grain size [µm]	30	[10 ⁻⁶ K ⁻¹]	14,2 (14,0)
Vickers hardness (HV 5)	* 205/225/230	Heat treatment	600 °C 15 min.
Elongation limit (Rp 0.2) [Mpa]	* 460/510/530	Soft annealing	750 °C 10 min.
Ductile yield (A5) [%]	* 21/15/13	(then quenching in water at 20 °C)	
Modulus of elasticity approx. [MPa]	113,000	* soft / after firing / hardened	

Processing instructions

Modelling:

- Minimum metal thickness (after grinding): for ceramic veneering 0,4 mm, for acrylic veneering with retention pearls 0,3 mm.
- Connecting parts between the pontics should be as thick and high as possible (at least 3,5 mm high and 2,5 mm wide).
- Spruing of single crowns: provide casting reservoir.

Investing: Use phosphate-bonded crown and bridge investment materials (e. g. BellaStar, Bellavest®).

Casting and finishing: General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only.

- In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material.
- Use ceramic or graphite crucible. Ceramic crucible: sprinkle some Auromelt HF melting powder on the casting ingots.
- Continue to heat after complete melting of the casting pieces:

Flame melting:	5–8 seconds
HF induction heating:	approx. 3 seconds
Resistance heating:	approx. 2–3 minutes.

Use fine carbide, ceramically bonded stones or BEGO sintered diamond milling tools for finishing.

Ceramic: Only use veneering ceramics which are in accordance with DIN EN ISO 9693 with firing temperatures of up to approx. 980 °C and which are resistant to green discoloration by means of silver oxides! Always follow the ceramic manufacturer's instructions!

- Always blast the surface to be veneered (Korox® 110, 2 – max. 3 bar) and clean the frame thoroughly (steam clean or boil in Aqua dest.). Allow to cool down normally after firing.

Oxide firing:

- Maintain oxide firing at 960 °C without vacuum for 10 minutes.
- The oxide can be blasted again prior to application of the ceramics.

Soldering:

- Support object in a soldering block of Bellatherm®. Prepare a gap of max. 0.2 mm with parallel walls.
- Soldering before firing with the flame (1125 °C): BegoStar®-Solder (Order no. 61081) and Minoxid or Fluxsol.
- Soldering after firing in furnace (810 °C): BEGO-Gold-Solder I (Order no. 61017) and Minoxid. Allow to cool normally.
- Acid-treat residual flux and metal oxides in Aurocid (60 °C, 20–30 seconds). Clean thoroughly (steam clean or boil in aqua dest.).

Laser welding: Filler material: Bio PontoStar®-wire diameter 0,4 mm (Order no. 61157) or accurately fitting, custom-ground casting pieces.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products undergo constant further development and are therefore subject to modification regarding design and composition.

Secondary effects: Like e. g. allergies to contents of the alloy or electrochemically based reactions are may very rarely occur.

Reciprocal actions: In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known incompatibilities and allergies to contents of the alloy.