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The enclosed BEGO-GOLD card is part of these processing instructions!

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Alloy characteristics

Biocertificate	<input checked="" type="checkbox"/>	Modulus of elasticity approx. [MPa]	90,000
Type	2	Melting interval [°C]	1005–935
BEGO-GOLD colour code	yellow 2	Casting temperature [°C]	1200
Density [g/cm ³]	16.5	Preheating temperature [°C]	700
Average grain size [µm]	40	Heat treatment	400 °C 15 min.
Vickers hardness (HV 5)	* 120/150	Soft annealing	750 °C 10 min.
Elongation limit (Rp 0.2) [MPa]	* 200/245	(then quenching in water at 20 °C)	
Ductile yield (A5) [%]	* 26/12	* Soft / after casting	

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
**Processing instructions****Modelling:**

- Minimum metal thickness (after grinding): for acrylic veneering with retention pearls 0.3 mm.
- Spruing of single crowns: provide casting reservoir.

Investing: Use special inlay investment material Auroplus® or plaster-bonded investment materials.

Casting and finishing: General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only.

- In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material.
- Use ceramic or graphite crucible. Ceramic crucible: sprinkle some Auromelt HF melting powder on the casting ingots.
- Continue to heat after complete melting of the casting pieces: Flame melting, immediate casting! HF induction heating, immediate casting! Resistance heating 1–2 minutes.

If Auroplus® is mixed with water, residual investment material can be brushed off easily under running water after deflasking. Then strip or blast with Korox® 25, Korox® 50 or Korox® 110 at 2–3 bar or blast-polish with Perlablast® micro or Perlablast® 125. If required, carbide, ceramically bonded stones or BEGO sintered diamond milling tools can be used for finishing.

Acrylic veneering: To simplify rubber polishing, blast-polish with Perlablast® micro (lead-free soda glass). Then rubber polish with BEGO's green rubber polishers for precious metal and polish with BEGO paste for preliminary and subsequent polishing of precious metals.

The relevant instructions must be followed for processing the acrylic veneering materials.

Soldering and polishing:

- Support object in a soldering block of Bellatherm®. Prepare a gap of max. 0.2 mm with parallel walls.
- Soldering with: BEGO-Gold Solder I (810 °C) / Minoxyd flux, BEGO-Gold Solder II (765 °C) / Minoxyd flux or Pre-flux U Gold Solder (810 °C) with pressed-in flux.
- Acid-treat residual flux and metal oxides in Aurocid (60 °C, 20–30 seconds). Clean thoroughly (steam clean or boil in aqua dest.).
- Polishing is possible in commercially available gold polishing units.

Laser welding: Filling material: InLloyd® 100-wire, dia. 0.4 mm (Order no. 61163) or cast one yourself (Acrylic profile dia. 0.45 mm, Order no. 53005).

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products undergo constant further development and are therefore subject to modification regarding design and composition.

Secondary effects: Allergies to contents of the alloy or electrochemically based reactions are may very rarely occur .

Reciprocal actions: In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known allergies to contents of the alloy.