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
**Modelling/Sprue system:** Minimum metal thickness (after grinding): for ceramic veneering 0.3 mm, for acrylic veneering with retention pearls 0.3 mm. Shape groove in cervical and palatal area. Avoid sharp edges. In the case of work on long bridges, divide plastic hollow stick between the front teeth and in the area around the canine teeth.

**Investing/preheating:** Use phosphate-bonded crown and bridges investment materials (e. g. Bellavest®). Follow processing instructions! Preheating temperature 900–1000 °C.

**Melting/casting:** General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each charge cast new metal only. In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material. Use only ceramic crucibles.

**Moment for casting:** Vacuum pressure casting with induction heating (Nautilus®) and centrifugal casting with induction heating (Fornax®). When the last solid component has submerged completely in the melt, continue to heat for 0 to 12 seconds depending on the induction capacity of the casting unit, then trigger casting. Follow operating instructions for Fornax® and Nautilus®. Flame centrifugal casting (Fundor): Cast when the last solid component has submerged completely in the melt and the melt clearly moves due to the flame pressure.

**After deflasking:** Blast residual investment material with Korox® 250 at 3–4 bar. Use fine carbide, ceramically bonded stones or BEGO sintered diamond milling tools for finishing.

**Ceramic:** Use ceramics according to ISO 9693 with firing temperatures of up to approx. 980 °C (e. g. VMK 95). Also suitable: Ceramics with reduced firing temperature (e. g. Omega 900, VM 13). Always follow the ceramic manufacturer's instructions! Always blast the surface to be veneered (pure Korox® 250: 3–4 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.). After cleaning, hold with artery clamps and refrain from touching.

**Oxide firing:** Not necessary. If oxide firing is carried out to check the surface (950–980 °C/5 minutes with vacuum), always blast the oxide again (with pure Korox® 250: 3–4 bar). Clean thoroughly (steam clean or boil in aqua dest.).

**Ceramic firing:** Always apply basic material in two firing operations. The first coating thin (washbrand), the second coating opaque. Wash off frame under running water before application of next ceramic coating. Long-term cooling recommended (cooling phase up to approx. 600 °C). Remove ceramics only mechanically. Hydrofluoric acid (HF) corrodes the metal frame.

**Acrylic veneering:** The corresponding instructions of the manufacturers must be followed when processing the veneering systems.

**Final work:** Blast visible metal surfaces with Korox® 50, then blast-polish external surfaces with Perlablast®. Rubber-polish with BEGO rubber polishers, and after that finish-polish with BEGO cobalt-chrome polishing paste (blue). Clean thoroughly (steam clean or boil in aqua dest.).

**Soldering:** Soldering prior to firing with the flame: Wirobond® solder (REF 52622) and Fluxsol flux (REF 52531). Soldering after firing in the furnace: WGL solder (REF 61079) and Minoxid flux (REF 52530). Long-term cooling recommended (cooling phase up to approx. 600 °C).

**Laser welding:** Filler material: Wiroweld wire Ø 0.35 mm (REF 50003) or Ø 0,5 mm (REF 50005).

**Secondary effects:** Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

**Reciprocal actions:** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Reactions:** In case of known incompatibilities and allergies to contents of the alloy.

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.