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Melting interval: 1100–1150 °C
Working temperature: approx. 1180 °C

Flux: Fluxsol

Soldering should be carried out with an apparatus which produces a precisely adjustable point flame, such as the Multiplex micro-flame soldering apparatus. If there is no soldering block necessary the object should be clamped into a soldering rest.

When soldering in a soldering block keep the gap about 0.2 mm wide. The points to be soldered must be exposed so that the flame can reach them satisfactorily. Make the soldering block up out of the soldering investment material Bellatherm® and water and pre-dry it in a furnace at about 300 °C. Do not use soldering investment containing plaster.

When soldering to repair defects, enlarge the defect to make it funnel shaped and apply Fluxsol to it. Break off an appropriately sized piece of Wirobond®-solder and coat it with Fluxsol. Place Wirobond®-solder in the place to be soldered and heat by moving the flame around over it. The solder can also be used for spot welding (spot welding apparatus). Before soldering coat with Fluxsol. As soon as the object is red hot direct the point flame onto the solder until it fills up the place to be soldered.

The Wirobond® soldering strips which are adapted to the casting material have the comparable characteristics as the BEGO-Co-Cr metal-to-ceramic alloys (Wirobond®), so that ceramic bonds equally well on the soldered places.

Secondary effects: Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

Reciprocal actions: In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Reactions: In case of known incompatibilities and allergies to contents of the alloy. Improper handling may pose a health risk.

Warranty: Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.

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