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**Modelling/Sprue system:** Always place sprues in the most solid wax-up areas, e. g. at the transition between saddle and base. Provide solid places which the melt can only reach through a thinly modelled area with an additional sprue (Ø 3 mm).

**Investing/preheating:** Use phosphate-bonded partial-denture investment materials (Wirovest®, WiroFine). Follow processing instructions!  
Preheating temperature 950–1050 °C.

**Melting/casting:** General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only. In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material. Use only ceramic crucibles.

**Moment for casting:** Vacuum pressure casting with induction heating (Nautilus®) and centrifugal casting with induction heating (Fornax®). When the last solid component has submerged completely in the melt, continue to heat for 2 to 12 seconds depending on the induction capacity of the casting unit, then trigger casting. Follow operating instructions for Fornax® and Nautilus®. Flame centrifugal casting (Fundor): Cast when the last solid component has submerged completely in the melt and the melt clearly moves due to the flame pressure.

**After deflasking:** Blast with Korox® 250 at approx. 4 bar. Critical areas – e. g. inner clasp sides and stress breakers – are to be blasted extremely carefully (Duostar or EasyBlast blasting units, Korox® 50 blasting material). Use fine carbide, ceramically bonded stones or BEGO sintered diamond milling tools for finishing. Polishing (Eltropol polishing unit, Wirolyt polishing liquid), rubber-polishing (BEGO rubber polisher, black) and finish-polishing (BEGO cobalt chrome polishing paste, blue). Clean thoroughly (steam clean or boil in aqua dest.).

**Soldering:** BEGO cobalt chrome soldering rods (REF 52520) and Fluxsol flux (REF 52531) are recommended. If soldering is carried out in combination with precious metal: Minoxid flux (REF 52530).

**Laser welding:** Filler material: Wiroweld wire Ø 0.35 mm (REF 50003) or Ø 0.5 mm (REF 50005).

**Secondary effects:** Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

**Reciprocal actions:** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Reactions:** In case of known incompatibilities and allergies to contents of the alloy.

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.