












BEGO compatibility overview 3D printing system components

Select your printer manufacturer

| | | | |
|---|---|---|---|
|  |  |  |  |
|  |  |  |  |
|  |  |  | <p>DLP printers in general</p> |

BEGO compatibility overview 3D printing system components



BEGO Varseo XS

Cleaning

| | | | | | |
|--|--|---|--|---|---|
|  VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix ¹ Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar ¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape ¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| |  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %) | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| |  VarseoSmile Teeth Denture teeth | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %) | | |
| |  VarseoWax CAD/CAST Burnout objects | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | |
|  VarseoWax Model Dental models | Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 5 min (Isopropanol 99 %) | Anycubic ¹ Wash & Cure Plus ³: 8 min (Isopropanol 99 % or Ethanol 96 %) | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax Surgical Guide Surgical guides and place- ment aids for implant prosthetics | | | | | |
|  VarseoWax Tray Individual impression trays | | | | | |

NEW!
Also available in
Bleach

Light-curing

| | | | | | | |
|--|---|--|---|---|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | SprintRay ProCure: 2 × 20 min @ 20 °C | | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A | Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n | |
| BEGO Otoflash: 2 × 500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Ackuretta CURIE/ Varseo Cure²: 1 × 2 min Exposure para- meters: P13 D8 T2.00 B0n | | | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| Flash or LED light-curing device, e.g.: | | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer): 2 × 180 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | Anycubic Wash & Cure Plus³: 2 × 20 min | Ackuretta CURIE/Varseo Cure²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹² PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



BEGO Varseo/Varseo L/Varseo S

Cleaning

| | | | | | |
|---|--|--|---|---|---|
|  VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix ¹ Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar ¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape ¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Teeth Denture teeth | | | | | |
|  VarseoWax CAD/CAST Burnout objects | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax Model Dental models | | | | | |
|  VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | Ultrasonic bath: 3 min + 2 min (Ethanol 96 %) | | | | |
|  VarseoWax Tray Individual impression trays | Ultrasonic bath: 3 min + 2 min (Ethanol 96 %) | | | | |

Light-curing

| | | | | | | |
|---|--|--|---|--|--|--|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox ¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| | | | | | | |
| BEGO Otoflash: 2 × 500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Ackuretta CURIE / Varseo Cure ²: 1 × 2 min Exposure parameters: P13 D8 T2.00 B0n | | | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| BEGO Otoflash: 2 × 1.000 + 2 × 2.000 flashes | HiLite Power (Kulzer ¹): 1 × 90 + 2 × 180 sec | | | | | |
| | | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer ¹): 2 × 180 sec | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹² PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.


BEGO compatibility overview 3D printing system components



Ackuretta ¹ DENTIQ / Freeshape 120 ⁴

Cleaning

Light-curing

| | | | | | | | | | | | | |
|--|---|---|--|--|--|--|---|--|--|--|--|--|
|  <p>VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | <p>Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %)</p> | <p>Whip Mix ¹ Veriwash + Veriwhirl / Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %)</p> | <p>Manual cleaning with tensides: See last page</p> | <p>SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer ¹): 2 × 90 sec</p> | <p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p> | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A</p> | <p>CUREbox ¹ Plus: 2 × 20 min @ 30 °C</p> | <p>Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | <p>Formlabs Form Wash: 3 min (Isopropanol 99 %)</p> | <p>Whip Mix Veriwash + Veriwhirl / Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %)</p> | <p>Manual cleaning with tensides: See last page</p> | <p>SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer ¹): 2 × 90 sec</p> | <p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p> | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p> | <p>CUREbox Plus: 2 × 20 min @ 30 °C</p> | <p>Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoSmile Teeth Denture teeth</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | | <p>Whip Mix Veriwash + Veriwhirl / Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %)</p> | | <p>SprintRay ProWash / Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer ¹): 2 × 90 sec</p> | | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A</p> | | <p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoWax CAD/CAST Burnout objects</p> | | | | | | | | | | | | |
|  <p>VarseoWax Model Dental models</p> | <p>Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %)</p> | <p>Formlabs Form Wash: 5 min (Isopropanol 99 %)</p> | <p>Anycubic ¹ Wash & Cure Plus ³: 8 min (Isopropanol 99 % or Ethanol 96 %)</p> | | <p>Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 2.000 flashes</p> | <p>HiLite Power (Kulzer ¹): 2 × 180 sec</p> | <p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p> | <p>Anycubic Wash & Cure Plus ³: 2 × 20 min</p> | <p>Ackuretta CURIE / Varseo Cure ²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n</p> | <p>Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum)</p> | |
|  <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p> | | | | | | | | | | | | |
|  <p>VarseoWax Tray Individual impression trays</p> | | | | | | | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
⁴ Ackuretta Dentiq & Freeshape 120 can only be used with anodized aluminium build platform (for VarseoSmile Crown ^{plus} and VarseoSmile Temp).
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹² PrograPrint object holder must be purchased separately (not supplied as standard).
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.


BEGO compatibility overview 3D printing system components



Ackuretta¹ SOL⁵

Cleaning

Light-curing

| | | | | | | | | | | | | |
|---|---|--|---|--|---|--|---|--|---|---|--|---|
|  <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | <p>Formlabs¹ Form Wash: 3 min (Isopropanol 99 %)</p> | <p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99 %)</p> | <p>Manual cleaning with tensides: See last page</p> | <p>SprintRay¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer¹): 2 × 90 sec</p> | <p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p> | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p> | <p>CUREbox¹ Plus: 2 × 20 min @ 30 °C</p> | <p>Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | <p>Formlabs Form Wash: 3 min (Isopropanol 99 %)</p> | <p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99 %)</p> | <p>Manual cleaning with tensides: See last page</p> | <p>SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer¹): 2 × 90 sec</p> | <p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p> | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p> | <p>CUREbox Plus: 2 × 20 min @ 30 °C</p> | <p>Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoSmile Teeth Denture teeth</p> | <p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %)</p> | | <p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99 %)</p> | | <p>SprintRay ProWash / Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99 %)</p> | <p>BEGO Otoflash: 2 × 1.500 flashes</p> | <p>HiLite Power (Kulzer¹): 2 × 90 sec</p> | | <p>SprintRay ProCure: 2 × 20 min @ 20 °C</p> | <p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A</p> | | <p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p> |
|  <p>VarseoWax CAD/CAST Burnout objects</p> | | | | | | | | | | | | |
|  <p>VarseoWax Model Dental models</p> | | | | | | | | | | | | |
|  <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p> | | | | | | | | | | | | |
|  <p>VarseoWax Tray Individual impression trays</p> | | | | | | | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

⁵ Ackuretta SOL can only be used with anodized aluminium build platform (Small, Medium and Large).

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹² PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.




BEGO compatibility overview 3D printing system components


Anycubic ¹ Photon Mono X ⁶

Cleaning

Light-curing



VarseoSmile Crown ^{plus}
Permanent single crowns,
inlays, onlays, and veneers



VarseoSmile Temp
Temporary crown and
bridge restorations, inlays,
onlays, and veneers



VarseoSmile Teeth
Denture teeth



VarseoWax CAD/CAST
Burnout objects



VarseoWax Model
Dental models



VarseoWax Surgical Guide
Surgical guides and place-
ment aids for implant
prosthetics



VarseoWax Tray
Individual impression trays

| | | | |
|---|--|--|---|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 5 min (Isopropanol 99 %) | Anycubic Wash & Cure Plus ³: 8 min (Isopropanol 99 % or Ethanol 96 %) | Rapid Shape ¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | |
| | | | |

| | | | | | |
|--|---|---|--|--|--|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| Flash or LED light-curing device, e.g.: | | | | | |
| BEGO Otoflash: 2 x 2.000 flashes | HiLite Power (Kulzer ¹): 2 x 180 sec | Formlabs Form Cure: 2 x 20 min @ 60 °C | Anycubic Wash & Cure Plus ³: 2 x 20 min | Ackuretta ¹ CURIE/Varseo Cure ²: 1 x 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) |
| | | | | | |
| | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
⁶ Compatibility applies to the design status to the serial no. P02123C0306508.
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Asiga¹ Max UV / Max 405

Cleaning

| | | | | | |
|--|--|---|--|--|--|
|  VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Teeth Denture teeth | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | Whip Mix Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | | SprintRay ProWash / Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax CAD/CAST Burnout objects | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax Model Dental models | Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 5 min (Isopropanol 99 %) | Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99 % or Ethanol 96 %) | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | | | | | |
|  VarseoWax Tray Individual impression trays | | | | | |

NEW!
Also available in Bleach

Light-curing

| | | | | | | |
|--|---|--|---|---|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 × 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 × 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A | | Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| BEGO Otoflash: 2 × 500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Ackuretta CURIE / Varseo Cure²: 1 × 2 min Exposure parameters: P13 D8 T2.00 B0n | | | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| Flash or LED light-curing device, e.g.: | | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer): 2 × 180 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | Anycubic Wash & Cure Plus³: 2 × 20 min | Ackuretta CURIE / Varseo Cure²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| | | | | | | |


¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹² PrograPrint object holder must be purchased separately (not supplied as standard).
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components




Cleaning

Light-curing



VarseoSmile Crown ^{plus}
Permanent single crowns,
inlays, onlays, and veneers



VarseoSmile Temp
Temporary crown and
bridge restorations, inlays,
onlays, and veneers



VarseoSmile Teeth
Denture teeth



VarseoWax CAD/CAST
Burnout objects



VarseoWax Model
Dental models



VarseoWax Surgical Guide
Surgical guides and place-
ment aids for implant
prosthetics



VarseoWax Tray
Individual impression trays

| | | | |
|---|---|---|--|
| Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 5 min (Isopropanol 99 %) | Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99 % or Ethanol 96 %) | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|---|---|---|--|

| | | | | | |
|--|--|---|---|--|--|
| Flash or LED light-curing device, e.g.: | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer¹): 2 × 180 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | Anycubic Wash & Cure Plus³: 2 × 20 min | Ackuretta¹ CURIE/Varseo Cure²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) |




¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
⁷ Printer must be operated in 4K mode.
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.



BEGO compatibility overview 3D printing system components

DEKEMA ¹ trix print ²

Cleaning

| | | | | | |
|---|---|--|---|---|---|
|  VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix ¹ Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar ¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape ¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | | |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | | |
|  VarseoSmile Teeth Denture teeth | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99 %) | | SprintRay ProWash / Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | | |
| | | | | | |

NEW!
Also available in
Bleach

Light-curing

| | | | | | | |
|---|--|---|--|---|--|--|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox ¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 × 2,5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A | | Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | | | |
| | | | | | | |
| | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹² PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.


BEGO compatibility overview 3D printing system components

Formlabs¹ Form 2⁸ (stainless steel build platform necessary)



Cleaning


Light-curing



VarseoSmile Crown^{plus}
Permanent single crowns, inlays, onlays, and veneers

| |
|--|
| |
|--|

| |
|--|
| |
|--|



VarseoSmile Temp
Temporary crown and bridge restorations, inlays, onlays, and veneers

| |
|--|
| Formlabs Form Wash: 3 min (Isopropanol 99 %) |
|--|

| |
|---|
| Formlabs Form Cure: 2 × 20 min @ 60 °C |
| Formlabs Fast Cure¹⁰: 2 × 2.5 min with intensity lvl 1 (see preprogrammed resin profile Temporary CB) |



VarseoSmile Teeth
Denture teeth

| |
|--|
| |
|--|

| |
|--|
| |
|--|



VarseoWax CAD/CAST
Burnout objects

| |
|--|
| |
|--|

| |
|--|
| |
|--|



VarseoWax Model
Dental models

| |
|--|
| |
|--|

| |
|--|
| |
|--|



VarseoWax Surgical Guide
Surgical guides and placement aids for implant prosthetics

| |
|--|
| |
|--|

| |
|--|
| |
|--|



VarseoWax Tray
Individual impression trays

| |
|--|
| |
|--|

| |
|--|
| |
|--|

Formlabs Form 3B/3B+ (stainless steel build platform necessary)



Cleaning

Light-curing

| |
|--|
| Formlabs Form Wash: 3 min (Isopropanol 99 %) |
|--|

| |
|--|
| Formlabs Form Cure: 2 × 20 min @ 60 °C |
| Formlabs Fast Cure¹⁰: 2 × 2.5 min with intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |



| |
|--|
| Formlabs Form Wash: 3 min (Isopropanol 99 %) |
|--|

| |
|---|
| Formlabs Form Cure: 2 × 20 min @ 60 °C |
| Formlabs Fast Cure¹⁰: 2 × 2.5 min with intensity lvl 1 (see preprogrammed resin profile Temporary CB) |

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
⁹ VarseoSmile Crown^{plus} is distributed by Formlabs as Permanent Crown + VarseoSmile Temp as Temporary CB.
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹³ Bleach shade only available for Form 3B/+
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Ivoclar PrograPrint PR5

Cleaning

| | | | | | |
|---|--|--|---|--|--|
|  VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |

NEW!
Also available in Bleach

Ivoclar ¹ PrograPrint Clean:
3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)

Rapid Shape ¹ RS Wash:
4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)



Ivoclar PrograPrint Clean:
3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %)

Rapid Shape RS Wash:
4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %)

Light-curing

| | | | | | | |
|--|--|--|---|--|--|--|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox ¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |

| | | | | | | |
|--|--|--|---|---|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer ¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure ²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure ¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |



¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹² PrograPrint object holder must be purchased separately (not supplied as standard).
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Microlay¹ Versus 385

Cleaning

| | | | | | |
|--|--|---|--|--|--|
|  VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Teeth Denture teeth | | | | | |
|  VarseoWax CAD/CAST Burnout objects | | | | | |
|  VarseoWax Model Dental models | Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 5 min (Isopropanol 99 %) | Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99 % or Ethanol 96 %) | | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | | | | | |
|  VarseoWax Tray Individual impression trays | | | | | |

Light-curing

| | | | | | | |
|--|--|--|---|--|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| | | | | | | |
| | | | | | | |
| Flash or LED light-curing device, e.g.: | | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer¹): 2 × 180 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | Anycubic Wash & Cure Plus³: 2 × 20 min | Ackuretta CURIE/Varseo Cure²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹² PrograPrint object holder must be purchased separately (not supplied as standard).
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.



BEGO compatibility overview 3D printing system components

Rapid Shape D10+ / D20 II / D20+ / D30 II / D30+

Cleaning

| | | | | | | | |
|---|---|--|---|---|---|--|--|
| VarseoSmile Crown ^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs ¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta ¹ Cleani ²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay ¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | |
| | | | | Ivoclar ¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape ¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | |
| VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani ²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | |
| VarseoSmile Teeth Denture teeth | | | | | | | |
| VarseoWax CAD/CAST Burnout objects | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | | | | | |
| | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | | | | |
| VarseoWax Model Dental models | Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 5 min (Isopropanol 99 %) | Anycubic ¹ Wash & Cure Plus ³: 8 min (Isopropanol 99 % or Ethanol 96 %) | | | | |
| | | | | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | | | |
| VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | | | | | | | |
| VarseoWax Tray Individual impression trays | | | | | | | |

NEW!
Also available in Bleach

Light-curing

| | | | | | | |
|--|---|--|---|---|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer)¹: 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 x 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| | | | | | | |
| BEGO Otoflash: 2 × 500 flashes | HiLite Power (Kulzer): 2 × 90 sec | Ackuretta CURIE/ Varseo Cure²: 1 × 2 min Exposure para- meters: P13 D8 T2.00 B0n | | | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| Flash or LED light-curing device, e.g.: | | | | | | |
| BEGO Otoflash: 2 × 2.000 flashes | HiLite Power (Kulzer): 2 × 180 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | Anycubic Wash & Cure Plus³: 2 × 20 min | Ackuretta CURIE/Varseo Cure²: 1 × 5 min Exposure parameters: P13 D8 T5.00 B0n | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | |
| | | | | | | |
| | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹² PrograPrint object holder must be purchased separately (not supplied as standard).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



SprintRay¹ Pro 95/Pro 95 S/Pro 55/Pro 55 S

Cleaning

| | | | | | |
|--|--|---|--|--|--|
|  VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar¹ PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
| | | | | Ivoclar PrograPrint Clean: 3 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + brush cleaning + spraying off (Isopropanol 99 %) | Rapid Shape RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoSmile Teeth Denture teeth | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | | Whip Mix Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99 %) | | SprintRay ProWash / Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99 %) |
|  VarseoWax CAD/CAST Burnout objects | | | | | |
|  VarseoWax Model Dental models | | | | | |
|  VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | | | | | |
|  VarseoWax Tray Individual impression trays | | | | | |

Light-curing

| | | | | | | |
|---|---|---|--|---|--|---|
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 × 2.5 min at intensity lvl 1 (see preprogrammed resin profile Permanent Crown) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE/ Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | Ivoclar PrograPrint Cure: 2 × 2 min on the PrograPrint object holder ¹² (25 % intensity, 405 nm) | Rapid Shape RS Cure: 15 min using lower and upper wavelength at 100 % Power (without vacuum) | Formlabs Fast Cure¹⁰: 2 × 2.5 min at intensity lvl 1 (see preprogrammed resin profile Temporary CB) |
| BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Teeth (60 sec – pause – 50 sec), Zone A | | Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |
| | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.
¹² PrograPrint object holder must be purchased separately (not supplied as standard).
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Whip Mix¹ VeriBuild / VeriEko⁹

Cleaning

Light-curing

| | | | | | | | | | | | | |
|--|--|---|--|--|---|--|---|--|---|--|--|--|
|  VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl / Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay¹ ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A | CUREbox¹ Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
|  VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers | Ultrasonic bath: 3 min + 2 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs Form Wash: 3 min (Isopropanol 99 %) | Whip Mix Veriwash + Veriwhirl / Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%) | Manual cleaning with tensides: See last page | SprintRay ProWash / Dry: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) | BEGO Otoflash: 2 × 1.500 flashes | HiLite Power (Kulzer¹): 2 × 90 sec | Formlabs Form Cure: 2 × 20 min @ 60 °C | SprintRay ProCure: 2 × 20 min @ 20 °C | SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A | CUREbox Plus: 2 × 20 min @ 30 °C | Ackuretta CURIE / Varseo Cure²: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n |
|  VarseoSmile Teeth Denture teeth | | | | | | | | | | | | |
|  VarseoWax CAD/CAST Burnout objects | | | | | | | | | | | | |
|  VarseoWax Model Dental models | | | | | | | | | | | | |
|  VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics | | | | | | | | | | | | |
|  VarseoWax Tray Individual impression trays | | | | | | | | | | | | |

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

¹⁰ Device needs to cool down between postcuring cycles for at least 10 minutes. This can be accelerated by blowing cold air into the device.

¹¹ Whip Mix VeriBuild & VeriEKO can only be used with anodized aluminium build platform (VeriEKO anodized aluminium build platform in Small, Medium and Large).

¹² PrograPrint object holder must be purchased separately (not supplied as standard).


For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components


DLP printers in general (with 385 – 405 nm wave length)

Cleaning

Light-curing



VarseoSmile Crown^{plus}
Permanent single crowns,
inlays, onlays, and veneers



VarseoSmile Temp
Temporary crown and
bridge restorations, inlays,
onlays, and veneers



VarseoSmile Teeth
Denture teeth



VarseoWax CAD/CAST



VarseoWax Model



VarseoWax Surgical Guide
Surgical guides and placement aids for implant
prosthetics



VarseoWax Tray
Individual impression trays

| | | | |
|---|---|---|--|
| Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %) | Formlabs¹ Form Wash: 5 min (Isopropanol 99 %) | Anycubic¹ Wash & Cure Plus²: 8 min (Isopropanol 99 % or Ethanol 96 %) | Rapid Shape¹ RS Wash: 4 min reservoir 1 (Isopropanol 99 %) + 3 min reservoir 2 (Isopropanol 99 %) + 3 min drying + spraying off (Isopropanol 99 %) |
|---|---|---|--|

Flash or LED light-curing device, e.g.:

| | |
|----------------------|-------------------------|
| BEGO Otofash: | HiLite Power |
| 2 × 2.000 | (Kulzer ¹): |
| flashes | 2 × 180 sec |

**Formlabs Form
Cure:**
2 × 20 min @
60 °C

Ackuretta¹ CURIE/Varseo Cure²:
1 × 5 min
Exposure parameters:
P13 D8 T5.00 B0n

Rapid Shape RS Cure:
15 min using lower and
upper wavelength at 100 %
Power (without vacuum)

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
Not all products and services shown are available in all countries.

3. Compatibility applies to the design status to the serial no. W31126A0405446.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

Manual cleaning with Tensides

This cleaning method is valid for VarseoSmile Crown^{plus} and VarseoSmile Temp

Necessary tools, equipment and materials

- InovaPrint wash (hp-dent¹) general purpose 3D print cleaner
- Tap water
- 1-propanol (70 Vol.-%)
- Toothbrush
- Brush
- Instrument for holding the printed objects (e.g. artery clamp)
- Absorbent pad (e.g. paper towels) or tub
- 2 cups for 1-propanol (70 Vol.-%)
- 1 cup for cleaning solution
- Compressed air with trigger/splash guard
- Personal protective equipment: protective gloves and goggles



Cleaning process

Carefully remove the printed object from the platform with the help of a spatula.



Clean the printed object by using a toothbrush and/or brush in two steps:

Step 1:

Pre-wash with reusable cleaning solution: 5 % InovaPrint wash (hp-dent) + 95 % tap water.

- 1.1 Vigorously swirl object in cleaning solution for 15 s using pliers or an arterial clamp. Ensure not to damage the printed object.



- 1.2 Remove excess resin by using a toothbrush and cleaning solution. In addition, a brush can be used to clean the inside of crowns. Shortly swirl object in cleaning solution as needed.



Note: Clean toothbrush and brush regularly with 1-Propanol when covered with resin. Dry before reusing toothbrush or brush again for cleaning the printed object.

- 1.3 Use compressed air under a fume hood with splash guard to remove the cleaning solution from the printed object's surface.



- 1.4 Repeat step 1.2 and 1.3 until just a thin layer of resin remains on the surface of the print.

Step 2:

Finish with fresh 1-Propanol (70 Vol.-%).

- 2.1 Vigorously swirl printed object in 1-Propanol for 5 s and dry object immediately with compressed air.



- 2.2 Check for remaining resin (shiny spots). If surface of printed object is matte, remove support structure as described in the following step (step 2.3). Shiny spots can be removed by quickly brushing the surface with a 1-Propanol soaked brush. Immediately dry object with compressed air.

- 2.3 Remove the support structure with the help of a cutting wheel or side cutters. Ensure not to deform the printed object.



- 2.4 Quickly clean occlusal surface with a 1-Propanol soaked brush to remove any excess resin. Immediately dry printed object with compressed air.



Note: Contact time of printed object with 1-Propanol needs to be reduced to a minimum to avoid the formation of white spots on the surface.

- 2.5 Post-cure the dental objects by using a validated post-curing unit.

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.