

## Wiron<sup>®</sup> 99

# Premium NiCr alloy for metal-to-ceramic work or composite veneering – beryllium-free

- Reliable metal-ceramic bond minimises the risk of splitting off at a later date
- Low hardness of 195 HV10 for easy, quick trimming and polishing to a high-lustre
- Problem-free processing in all induction casting machines easy casting point recognition
- High degree of intraoral comfort for patients due to high thermal stability
- Biocompatible and corrosion-resistant



## The superior properties

Wiron<sup>®</sup> 99 has been processed successfully for over 30 years. The composition guarantees reliable recognition of the casting point both with centrifugal casting machines and vacuumpressure casting machines, e.g. Nautilus<sup>®</sup> or Fornax<sup>®</sup>

## The advantages for you:

- Easy casting
- Reliable further processing

The hardness of Wiron $^{\odot}$  99 (195 HV10) is the same as that of a high-gold-content alloy.

## The advantage for you:

- Easy, quick trimming
- Polishing to a high-lustre

The Young's modulus – decisive for the load-bearing capacity of crown and bridge alloys – is about twice as high as that of precious metal.

#### The advantages for you:

- The load-bearing capacity of a bridge is twice as high with the same design
- Reliable protection against deformation by masticatory forces

The high thermal stability provides a high degree of safety against deformation when welding, soldering and firing the porcelain.

#### The advantage for you:

- No distortion
- Perfectly fitting frameworks

The low thermal conductivity of Wiron<sup>®</sup> 99 protects the pulp of the abutment teeth against unpleasant temperature stimuli in the patient's mouth.

#### The advantage for you:

• Provides increased intraoral comfort, ensuring complete patient satisfaction

## **Precise fit**

Optimum castings in terms of accuracy, fit and surface quality can be achieved using Bellavest<sup>®</sup> SH, a BEGO precision casting investment. Bellavest<sup>®</sup> SH is a shock-heat investment for crowns and bridges that also produces precise castings with standard preheating.

#### The advantage for you:

• Quick, economic working because time-consuming adjustments are no longer required

## Reliable porcelain veneering

The outstanding bond strength between  $\mathsf{Wiron}^{\circledast}$  99 and porcelain is the result of close collaboration with leading porcelain manufacturer

### The advantages for you:

- Reliable veneering with all commercially available porcelains with a corresponding CTE value
- No additional, expensive bonder required
- Durable, reliable bond minimises the risk of splitting off intraorally at a later date

## **Biocompatible and corrosion-resistant**

Chrome and molybdenum ensure the excellent corrosion behaviour of Wiron<sup>®</sup> 99.

The advantages for you:

- No cytotoxic potential, no skin irritation or allergic sensitisation.
- Guarantee of the biocompatibility
- High degree of safety for the dentist and patient (biocertificate available at www.bego.com)



Wiron® 99 is the result of over 30 years of clinical experience and systematic development

## Product details

Alloy characteristics	Standard values	
• Alloy type (ISO 22674)	3	
• Density	8.3 g/cm <sup>3</sup>	
Preheating temperature	900–1000 °C	
Solidus, liquidus temperature	1310, 1360 °C	
Casting temperature approx.	1450 °C	
• Young's modulus	170 GPa	
• Proof strength ( $R_{p0.2}$ )	335 MPa	
• Ultimate strength (R <sub>m</sub> )	655 MPa	
• Elongation after fracture (A <sub>5</sub> )	43 %	
• Vickers hardness (HV10)	195 HV10	
<ul> <li>Coefficient of thermal expansion (CTE) 25–500 °C. 10<sup>-6</sup> K<sup>-1</sup></li> </ul>	13.9	

#### Composition in % by mass

• Ni 65.6 · Cr 22.5 · Mo 9.5 · Si 1.0 · Ce · Mn · Nb				
Accessories and	Presentation	Content	REF	
<ul> <li>Wiron<sup>®</sup> 99</li> <li>Wiron<sup>®</sup> 99</li> </ul>	1 Pack 1 Pack	1000 g 250 g	50225 50226	
• Wiroweld NC, Ni-Cr laser wire, carbon-fre	1 Roll	5.5 m – 4 g	50006	
Wiron <sup>®</sup> soldering rods	1 Pack	4 g	52625	
<ul> <li>Bellavest<sup>®</sup> SH</li> <li>1 carton 4.8 kg</li> <li>1 carton 12.8 kg</li> </ul>	160 g bags 160 g bags	30 pieces 80 pieces	54247 54252	
<ul> <li>BegoSol<sup>®</sup> HE mixing liquid 1 bottle 1 canister</li> </ul>		1   5	51095 51096	
ISO 22674 · ISO 9693-				